Hemolex[®] LA

38342-27



Orion Diagnostica

Explanation of symbols used on labelling Erklärung der Symbole • Explication des symboles

	English	Deutsch	Français	
	For <i>in vitro</i> diagnostic use	<i>In-vitro</i> - Diagnostikum	Diagnostic <i>in vitro</i>	
	Catalogue number	Bestell- Nummer	Référence	
	Batch code	Lot- Nummer	Numéro de lot	
	Use by	Verwendbar bis	à utiliser jusqu'à	
 Store at 2...8°C	Storage temperature	Lagertemperatur	Conservation	
	Consult instructions for use	Siehe Gebrauchs- information	Consulter la notice d'utilisation	
	Contents	Inhalt	Contenu	
	Contains sodium azide	Enthält Natriumazid	Contient de l'azide de sodium	
	Latex Reagent (vial 1)	Latex Reagenz (Fläschchen 1)	Réactif Latex (flacon 1)	
	Control Latex (vial 2)	Kontroll-Latex (Fläschchen 2)	Latex-Contrôle (flacon 2)	
	Positive Control (vial 3)	Positive Kontrolle (Fläschchen 3)	Contrôle Positif (flacon 3)	
	Origin: porcine / swine	Ursprung: Schwein	Origine: porc	
	Origin: human	Ursprung: human	Origine: humain	

Significado de los símbolos • Spiegazione dei simboli
Vysvětlivky použitých symbolů
Förklaring av symbolerna • Symbolien selitykset

	Español	Italiano	Česky	Svenska	Suomi
	Para uso diagnóstico <i>in vitro</i>	Per uso diagnostico <i>in vitro</i>	Diagnostický zdravotnický prostředek <i>in vitro</i>	För <i>in vitro</i> diagnostik	<i>In vitro</i> -diagnostiikkaan
	Número de catálogo	Codice	Katalogové číslo	Artikelnummer	Tuotenumero
	Número de lote	Numero di lotto	Číslo šarže	Batch nr.	Eränumero
	Fecha de caducidad	Utilizzare entro	Použitelné do	Används före	Käytettävä ennen
	Temperatura de conservación	Temperatura di conservazione	Teplotní rozmezí	Förvarings-temperatur	Säilytyslämpötila
	Consultar las instrucciones de uso	Consultare le istruzioni d'uso	Viz návod k použití	Se bruksanvisningen	Lue käyttöohjeet
	Contenido	Contenuto	Obsah	Innehåll	Sisältö
	Contiene azida de sodio	Contiene sodio azide	Obsahuje azid sodný	Innehåller natriumazid	Sisältää natriumatsidia
	Reactivo Latex (1)	Reagente al Lattice (1)	Latexová reagentie (1)	Latex reagens (flaska 1)	Lateksireagenssi (pullo 1)
	Control Latex (2)	Reagente di Controllo (2)	Kontrolní latexová reagentie (2)	Kontroll latex (flaska 2)	Kontrollilateksi (pullo 2)
	Control Positivo (3)	Controllo Positivo (3)	Pozitivní kontrola (3)	Positiv kontroll (flaska 3)	Positiivinen kontrolli (pullo 3)
	Origen: cerdo	Origine: maiale	Původ: prasečí	Ursprung: svin	Alkuperä: sika
	Origen: humano	Origine: umana	Původ: lidský	Ursprung: human	Alkuperä: humaani

1**2****3****4**

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- 1 Harmoinen A, Karvonen A-L, Koivula T, Winter M. Hemolex, a new latex test for occult blood in faeces, presented at the Nordiske Kongres i Klinisk Kemi, 13–15 August 1986, Denmark.
- 2 Winawer SJ, Fletcher RH, Miller L, Godlee F, Stolar MH, Mulrow CD, Woolf SH, Glick SN, Ganiats TG, Bond JH, Rosen L, Zapka JG, Olsen SJ, Giardiello FM, Sisk JE, van Antwerp R, Brown-Davies C, Marciniak DA, Mayer RJ. Colorectal cancer screening: Clinical Guidelines and rationale; *Gastroenterology* 1997;112:594–642.
- 3 Shike M, Winawer SJ, Greenwald PH et al. Primary prevention of colorectal cancer; the WHO Collaborating Centre for the Prevention of Colorectal Cancer, *Bulletin of the World Health organization*, 68(3):377–388, 1990.
- 4 Parkin P, Whelan SL, Ferlay J, Raymond L, Young J (1997). *Cancer incidence in five continents*, VOL VII IARC Scientific publication no. 143, 1997.
- 5 Mandel JS, Bond JH, Church TR, Snower DC, Bradley GM, Schuman LM, Ederer F. Reducing mortality from colorectal cancer by screening for fecal occult blood. Minnesota Colon cancer Study. *N Engl J Med* 1993; 328:1365–1371. (Published erratum appears in *N J Engl J Med* 1993;329:672.)
- 6 Kronborg O, Fenger C, Olsen J, Jørgesen OD, Søndergaard O. Randomised study of screening for colorectal cancer with faecal-occult-blood test. *Lancet* 1996 Nov 30; 348 (9040):1467–71.

- 7 Hardcastle JD, Chamberlain JO, Robinson MH, Moss SM, Amar SS, Balfour TW, James PD, Mangham CM. Randomised controlled trial of faecal-occult-blood screening for colorectal cancer. *Lancet* 1996 Nov 30;348 (9040) 1472–7.
- 8 Mandel JS, Church TR, Ederer F, Bond JH. Colorectal Cancer Mortality: Effectiveness of Biennial Screening for Fecal Occult Blood *J Natl Cancer Inst* 1999 Mar 3;91(5):434–7.
- 9 Nakama H, Kamijo N, Fujimori K, Fattah AS, Zhang B. Relation between fecal sampling times and sensitivity and specificity of immunochemical fecal occult blood tests for colorectal cancer: a comparative study. *Dis Colon Rectum* 1997 Jul;40(7):781–4.
- 10 Data on file, Orion Diagnostica.
- 11 Favenc L, Kapel N, Meiller D, Chochillon C, Gobert JG. Detection of occult blood in stools: comparison of three guaiac tests and a latex agglutination test, *Ann Biol Clin* 50 (1992), 311–313.
- 12 Goberts G: Comparison of two immunochemical methods for faecal occult detection, *Clin Chem News*, 18(6), (1992).

Hemolex® is a registered trademark of Orion Diagnostica Oy.

Manufacturer • Hersteller • Fabricant • Fabricante
Produttore • Výrobce • Tillverkare • Valmistaja:

Orion Diagnostica

Orion Diagnostica Oy
P.O.Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794
www.oriondiagnostica.fi